

APPLICATION OVERVIEW:

This application can be used to relieve pressure, provide support, and increase circulation to the area.

ITEMS NEEDED:

2 pre-cut strips of StrengthTape Scissors



Avoid touching the ends of the tape.

Clean any oils and

lotions form area.




After applying, firmly rub the tape to activate adhesive.



Trim excess hair for better adhesion.

SHIN SPLINTS INSTRUCTIONS





Always anchor the first and last two inches with absoultely no stretch.

STEP 1:



Anchor the strip at the top of the foot. Apply the strip up the shin, directly over the point of pain. Anchor the end with zero stretch.

BODY POSITION:



Bend leg to a 45 degree angle with toes pointed forward.

STEP 2:



Anchor one of the half strips just *above* the point of pain, and apply the strip horizontaly across the shin with 75% stretch. Anchoring with no stretch.

PREPERATION:



Cut a full strip in half with rounded corners to make 2 short strips.

STEP 3:



Repeat with the 2nd half strip, but this time anchoring just *below* the point of pain. Firmly rub the application to activate adhesive.

Please seek care if you have extreme pain, swelling, difficulty walking, or unnatural joint movements. These instructions are for educational purposes only and should not be used in place of professional medical treatment. Copyright Endevr LLC. All rights Reserved.