

FULL KNEE SUPPORT

AVERAGE DAYS WEAR PRO

APPLICATION OVERVIEW

Knee pain can be caused by any number of issues. The kneecap, or patella, could be moving incorrectly. One or more meniscus may be torn, ruptured, or inflamed. There may be arthritis, plica, chrondomalacia, or any number of issues with the bones that constitute the knee. KT Tape helps treat any of these conditions by relieving pressure, providing stability, and improving muscle function in the knee.

WHAT YOU NEED

2 strips of KT TAPE 1 full 10" strip 1 full 10" strip cut in half

BEFORE YOU START

APPLY BEFORE ACTIVITY Apply one hour before beginning activity

> **CLEAN SKIN** Clean dirt, oils and lotions from area

ACTIVATE ADHESIVE After application rub tape vigorously to activate adhesive






WATCH THE VIDEO kttape.com/instructions/full-knee-support BODY POSITION

Bend the knee at a 90 degree angle.



STRIP ONE



ANCHOR: Anchor middle of a half strip of tape under the kneecap with 80% stretch.



APPLY: Lay the ends down without stretch. A second half strip can be added for additional support.




STRIP TWO

FINISH: Anchor a full strip 10 inches above the kneecap.



ANCHOR: Apply the tape around the outer edge of the kneecap with the tape with 25% stretch.

STRIP THREE



APPLY: Apply tape across the lower knee with 50% stretch.



FINISH: Lav the last two inches down without stretch.



ANCHOR: Repeat the steps **APPLY:** Apply tape across for the previous strip with a the lower knee with 50% second full strip, mirroring stretch.



FINISH: Lay the last two inches down without stretch.

CAUTION: If you have skin sensitivities, cancer, or are pregnant, consult your doctor before use. Discontinue use if skin becomes irritated or sore. KT TAPE* is not a replacement for professional medical care. Warranties and remedies limited to product replacement cost. READ ALL CAUTIONS ON ENCLOSED INSTRUCTION SHEET PRIOR TO USE.

0% STRETCH

the application on the other side of the knee.