

WRIST

AVERAGE DAYS WEAR

APPLICATION OVERVIEW

Causes of wrist pain may include typing for long periods of time, overuse in sports, landing on an outstretched hand, or carrying heavy loads. KT Tape helps treat this condition by relieving pressure to reduce pain, relaxing associated muscles, and may increase circulation.

WHAT YOU NEED

2 strips of KT TAPE 1 full 10" strip 1 full 10" strip cut in half

BEFORE YOU START

APPLY BEFORE ACTIVITY Apply one hour before beginning activity

> CLEAN SKIN Clean dirt, oils and lotions from area

ACTIVATE ADHESIVE After application rub tape vigorously to activate adhesive






WATCH THE VIDEO kttape.com/instructions/wrist Hold arm out with hand pointing down

BODY POSITION





ANCHOR: Anchor the middle of a half strip of tape to the top of the wrist with 80% stretch.



FINISH: Lay ends down around the wrist without stretch.



ANCHOR: Anchor the middle of a second half strip to the bottom of the wrist as shown with 80% stretch.



FINISH: Apply the ends of tape without stretch around the wrist.

STRIP THREE



ANCHOR: Anchor a full strip on back of hand without stretch.



APPLY: Apply tape up the arm with 25% stretch.

FINISH: Lay the last two inches down without stretch.

CAUTION: If you have skin sensitivities, cancer, or are pregnant, consult your doctor before use. Discontinue use if skin becomes irritated or sore. KT TAPE* is not a replacement for professional medical care. Warranties and remedies limited to product replacement cost. READ ALL CAUTIONS ON ENCLOSED INSTRUCTION SHEET PRIOR TO USE.