

APPLICATION OVERVIEW:

This application can be used to relieve pressure, provide support, and increase circulation to the area.

ITEMS NEEDED:

2 pre-cut strips of StrengthTape Partner



Avoid touching the ends of the tape.

Clean any oils and

lotions form area.




After applying, firmly rub the tape to activate adhesive.



Trim excess hair for better adhesion.

HAMSTRING INSTRUCTIONS





Always anchor the first and last two inches with absoultely no stretch.

BODY POSITION:



Stretch the hamstring area as much as possible.

END RESULT:



Firmly rub to activate adhesive.

STEP 1:



Anchor the 1st strip above the inside back of the knee.

STEP 2:



Continue up the back of the thigh with a 25% stretch. Anchor off tape with zero stretch. STEP 3:



Mirror the 1st strip by anchoring the 2nd strip next to the 1st and continuing up the thigh with a 25% stretch.

Please seek care if you have extreme pain, swelling, difficulty walking, or unnatural joint movements. These instructions are for educational purposes only and should not be used in place of professional medical treatment. Copyright Endevr LLC. All rights Reserved.